

# Personal Protection Gown

## Protection, Comfort and Convenience for the Caregiver.

These gowns are the ideal impervious barrier protection for caregivers in the emergency room as well as in patient rooms. Thumbhooks prevent the gown from riding up and exposing wrists. The pullover design eliminates ties at the neck and the perforations in the back allow for fast removal. Gowns come in universal size and are made of latex-free, blue embossed polyethylene.

- Meets ASTM F1670-97 requirements for synthetic blood test penetration
- Meets ASTM F1671-97B requirements for viral penetration



### CATALOG NO. 235

Individually packed 15/box, 5 boxes/carton



#### THUMBHOOK STIRRUPS

Makes glove donning easier and protects against wrist exposure.



#### COMFORTABLE & CONVENIENT

Open back keeps wearer cool and comfortable. Perforations allow for fast, easy removal.



#### EASY TO USE PACKAGING

Individually packed 15/box. Hooks on to mounting bracket on most personal protection stations.

**busse**  
Hospital Disposables

P.O. Box 011067 Hauppauge, NY 11788-0920  
800.645.6526 • Fax: 631.435.4721  
www.busseinc.com • bussemail@aol.com

© 2000 Busse Hospital Disposables  
E.O.E. - M/F/D/V  
235